

DATASHEET

Name		PROGESTERONE	
Product	Progesterone	Article-No.	90133
Formulation	Off-white powder	CAS-No.	57-83-0
Mol.Weight	314.5	Formula	C ₂₁ H ₃₀ O ₂
Storage	Room temperature 15-30°C Protect from moisture		

- Conforms with BP2000 standard

Parameter	Specification
Assay	97 – 103 %
Specific optical rotation	186 - 194° (Ethanol)
Melting point	128 - 132°C
Related foreign steroids	< 1 %
Loss on drying	< 0.5 %