

## DATASHEET

<b>Name:</b>	<b>EGF, rec. Human</b>		
Product	Epidermal Growth Factor, rec. Human	Article-No.	40102
Formulation	solid, white	CAS-No.	62253-63-8
Storage	Freezer at -20°C Protect from sunlight		

Parameter	Specification
Purity	≥ 95 %
Activity	$3.5 \cdot 10^5 \sim 1 \cdot 10^6$ IU/mg (70 ~ 200% of the indicated amount $5 \cdot 10^5$ IU)
Identification	positiv
Solubility	1mg/1ml distilled H <sub>2</sub> O. It should not contain insoluble matter seen by naked eye
pH Value	7.2 ±0.2
Water residue	≤ 3.0 %
Sterility test	No bacterial growth